

UVCut Software manual v2.0

Molecular Machines & Industries AG (MMI)

www.molecular-machines.com

MMI CELLTOOLS**MMI UVCUT – User Manual:**

1	SECURITY ADVICE	5
1.1	Laser safety	5
1.2	General safety	5
2	INSTALLATION OF THE SYSTEM	6
3	THE MMI CELLTOOLS INSTRUMENTATION FAMILY	6
4	PRINCIPLES OF MMI CELLCUT®	7
4.1	System set up	7
5	HANDLING OF SAMPLES	8
5.1	Single Step collection using the MMI Isolation Cap®	8
5.2	Preparation of the slides	9
6	MMI CELLTOOLS® – SOFTWARE DESCRIPTION	10
6.1	Overview of functions	10
6.2	Quick Start	11
6.2.1	System start-up	11
6.2.2	System preparation	11
6.2.3	Handling a new slide	11
6.2.4	System turn off	12
6.3	MMI CELLTOOLS® – Main Application	13
6.3.1	Main Window and plugins	13
6.3.2	User specific database	13
6.3.3	Video settings	15
6.3.4	Camera configuration	15
6.3.5	Fluorescence	16
6.3.6	XY stage motion	17
6.3.7	Slide overview for navigation	17
6.3.8	Additional Settings for the xy-stage	19
6.3.9	Calibration procedures	19
6.3.9.1	Camera orientation	19
6.3.9.2	Objective calibration	20
6.3.9.3	Lens OffSet calibration	22
6.3.10	Multi-User Report	23

6.3.11	Hardware self monitoring	24
6.3.12	Overview Image Quality	24
6.4	MMI CELLTOOLS® – UV CUT® Plug-in	25
6.4.1	Drawing the cutting contour	25
6.4.2	Using the MMI IsolationCap®	27
6.4.3	Cutting of objects	28
6.4.4	Group handling	29
6.4.5	Distance measurement	30
6.4.6	Pin positions	30
6.4.7	Autodocumentation	30
6.5	Calibration	31
6.5.1	Laser position	31
6.5.2	Correction of slide offset	32
7	DISPOSABLES	33
7.1	Slides	33
7.2	Caps	33
7.3	Cell Chamber	33
8	MAINTENANCE	34
9	TROUBLE SHOOTING	35
9.1	View	35
9.2	Movement	35
9.3	Drawing	35
9.4	Cutting	35
9.5	Collection	36
10	WARRANTY	37
11	CUSTOMER SUPPORT	38
12	SYSTEM DATA OVERVIEW	39
A	SOFTWARE AND HOT KEY LIST	40
A.1	Shortcuts:	40

A.2 Hotkeys:**41**

1 Security advice

The system should only be used for microdissection as described in the manual. Do not use the system for any other purpose. Damage due to unauthorized use is not subject to warranties. Persons who have been properly trained should only use the system.

Read completely the Security Advice in Section 1 and the Manual before operation.

1.1 Laser safety

This system contains a laser for microdissection. The system includes safety devices to prevent laser interference with the user. Due to the power losses incurred as the beam passes through the microscope objectives the average power is nominal at the microscope work surface (i.e., at the membrane slide/tissue surface). If the illumination arm is and tilted backwards, the lasing action is interrupted by an electrical interlock. A blocking laser absorption filter is inside the microscope and tilted backwards, the lasing action is interrupted and provides the necessary eye protection for those times when the oculars may be used

To ensure full safety of the system please follow the steps below:

- Turn off the laser with the key switch to prevent unauthorized operation of the system.
- Do not remove any objective while the laser is operating.
- Never stare into the objective turret while the laser is operating.
- Use only the provided objectives for the laser microdissection.
- Never place reflecting objects in the beam path.
- The laser source and the optical equipment are enclosed within the blue housing. To avoid electrical or laser hazards, do not open the housing. By design, the laser beam is contained in a well-defined beam path, which is not serviceable, by the user.

1.2 General safety

Do not disassemble the system. The installation of the system is provided by MMI service personnel or MMI designated representative. Repairs, removal or exchange of components beyond the operations described in this manual may only be carried out by MMI service personnel or persons expressly authorized by MMI to do so. If you have any problems with the instrument, contact MMI.

The power supply is installed by MMI. MMI assures that the system is provided with the appropriate voltage. Do not change the power cords.

Avoid wet or dusty conditions near the system. If liquid gets inside the system, do not attempt to use it. Contact MMI.

Unplug all electrical supply before cleaning the system. Do not use cleaning fluids or sprays but only smooth and dry cloth.

If the stage control is not calibrated, table movements can be sudden and fast. Assure that the work area around the table is free of clutter and material.

Read the manual of your microscope for specific microscope precautions. If you do not have the manual contact your microscope provider or MMI.

2 Installation of the system

The system is installed by MMI Service personnel or our designated representative in the laboratory of the customer. After the installation training will be provided in the use and operation of the system. The customer should not change the installation of the equipment (see precautions in Section 1).

*Installation by
MMI*

3 The MMI CELLTOOLS instrumentation family

The MMI Cell Tools are a fully modular instrumentation family, including the following components:

- MMI CELLCUT® laser microdissection to isolate single cells or areas of tissue
- MMI CELLMANIPULATOR® optical tweezing to manipulate cells or beads with an optical trap
- MMI CELLECTOR® automated micro-pipetting to mechanically manipulate cells or beads with a capillary and mechanical micromanipulator
- MMI CELLEXPLORER pattern recognition software for PC based image analysis
- Fluorescence microscopy extension for both fluorescence and quantum dot applications

Any or all of these modules can be combined in one microscopic environment.

4 Principles of MMI CELLCUT®

The MMI CELLCUT® is a method to isolate, under microscopic view, small areas or single cells from histological sections for further microbiological analysis. Only the cell(s) wanted for further investigation are cut out. DNA, RNA, as well as proteins from undisturbed, pure samples can be investigated. No mechanical contact is necessary for the laser microdissection of the samples. Thus the method avoids contamination of the samples.

Isolation of cells

MMI CELLCUT® permits comfortable working with high precision and a large throughput without contamination.

Advantages of the system

Several areas of interest can be microdissected in one automated operation and collected in the MMI ISOLATIONCAP®. One MMI ISOLATIONCAP® used for Single Step Collection can collect several dissections, even from different slides.

MMI ISOLATIONCAP®

Review of the sample and selection of areas to cut do not have to be directly followed by the cutting operation. The automated cutting can be carried out later and by another person.


*Scan and cutting –
time shifted operation possible*

4.1 System set up

The MMI CELLCUT® system consists of a high performance research microscope with motorized scanning stage, an electronically controlled, solid-state laser, requisite laser beam delivery and transfer optics and a high-end Pentium computer with Windows and the sophisticated control software MMI CELLTOOLS.

Components

All usual microscope features are available. The MMI CELLTOOLS software controls the laser, image capture, and scanning stage actions without blocking any other microscope action.



In the figure,
change *micro-*
scopic objective to
microscope objec-
tive

5 Handling of samples

5.1 Single Step collection using the MMI Isolation Cap®

Microdissection for the isolation of cells is only useful when you can remove the parts of the tissue you are interested in from the tissue surroundings and from the slide. The single step collection makes sampling of one or several isolated areas easy and contamination free.

*Separation of cells
from the slide*

The Single Step Collection uses a protective membrane and a reaction tube with a special adherent lid.


*Single Step
Collection*

The purpose of the protective membrane is:

*Protective
membrane*

to avoid contamination of the sample and

to facilitate easy removal and collection of the cut area



With the laser beam originating from below the stage, the laser cuts through the tissue and the membrane. The separated tissue and membrane are collected with the lid of the reaction tube.

The laser cuts around the cells to be isolated. For the collection step, no additional radiation is used.

*No radiation
directly on
dissectates*

The membrane and the adhesive lid are chemically inert and have no influence on further molecular biological processing.

The membrane material is a thin PET membrane. The membrane is transparent and does not perturb the light beam.

The lid of the reaction tube also contains a diffuser insert. The lid improves the image quality remarkably and can be placed directly on top of tissue and membrane during all operations.

Image quality

5.2 Preparation of the slides

The Single Step Collection of the dissection requires special slides. The slides are provided with a 1.4µm thick, clear PET membrane.

Protective membrane


Intermediate slide



The tissue is mounted on the membrane as on an ordinary slide.

Paraffin sections, cryo sections or smears can be used.

Dissection slide



After the usual processing (deparaffinize staining, etc) the membrane with the tissue is inverted and placed onto a new glass slide and fixed in position on the microscope stage.

Thus the tissue is now under the membrane and protected against contamination.

6 MMI CELLTOOLS® – Software description

6.1 Overview of functions

The MMI CELLTOOLS® together with the MMI UVCUT® plug-in provides all necessary controls for:

- displaying live video
- saving images
- motorized xy-stage
- laser adjusting the power and laser focus
- storing the preferred camera settings
- manipulating the automatic cap lift

This results in an easy and user-friendly method for:

- scanning the sample
- documenting the sample
- marking the path for the laser cutting around single cells or cell clusters
- marking several areas to be cut in one operation
- storing the marked cutting paths (for later cutting)
- dissecting the marked sample areas automatically
- collecting the dissected areas without radiation radiation or risking contamination of the isolated material(s)

The operation of all modules is controlled by the comfortable and easy to use MMI CELLTOOLS®. The main application gives all the tools necessary for displaying the video, saving the image, adjusting the camera control and basic xy-stage movement.

The MMI UVCUT® plug-in adds the specific MMI CELLCUT® functions module.

Additional plugins for MMI CELLMANIPULATOR® and MMI CELLEXPLORER® are available.

*MMI CELLTOOLS®
software with
UVCUT® plug-in
Additional plugins*

6.2 Quick Start

6.2.1 System start-up

To start MMI CELLCUT®, proceed the following steps as follows:

- Start up the PC and allow the boot process to complete, reaching the Windows desktop.
- Turn on the microscope white light power supply.
- Energize the electronic controller with the key. A yellow LED will illuminate.
- Start MMI CELLTOOLS® software and wait until the software has finished the load procedure.
- Power the laser by pressing the button on the electronic box controller. A green LED will illuminate.

6.2.2 System preparation

For the handling of the MMI UVCUT®, start with the definition of a new setup.

For simplicity we call the first setup "Default". In this setup the laser power and focus settings, the on-screen laser position, the xy-stage calibration and the camera settings for each objective will be stored. To reach the setup explorer, using the PC mouse press "E" on the left side of the setup list.

If the objective you want to use does not appear in the objective list, define a new objective.

For each objective you prefer to use follow the calibration procedure for laser position, laser focus, laser power settings, xy-stage calibration (Chapter 6.3.9.2) and camera settings (Chapter 6.3.4).

6.2.3 Handling a new slide


Prepare your sample as described in Chapter 5 or following the detailed application notes provided by MMI.

Define a new slide by using the menu item:

File->New Slide

or press  on the left side of the setup list.

The slide editor will open and you can add a new slide by pressing the " + " button. All documentation will be stored in a




Define new setup



Choose objective



Prepare sample



Define new sample slide




folder with the same name, which you used for the slide.

Define the membrane slide work area by moving to the upper left corner of the membrane window and pressing  Limit1 and then by moving to the bottom right corner for the membrane window and pressing  Limit2.

Use the objective with the lowest magnification (standard is 4x). By pressing the scan button the software creates an overview of sample(e.g., a "roadmap").

By double clicking into the overview you can navigate into the area of interest. You can also use the cursor keys to move the xy-stage or drag the blinking rectangle with the mouse

Move the Caplift down by pressing the corresponding button. After lowering the Caplift you should adjust the microscope stage fine focus, as needed.

Activate any button in the toolbar with a mouse click. For example the freehand drawing tool  and enables the user to draw a line around the object to be microdissected.

Press the cut button to microdissect the marked object.

Lift the cap holder to collect the microdissected sample.

Now everything is clear to take out the cap and proceed with application.




Create overview



Navigate into the area of interest



Caplift down



Cut



Caplift Up

6.2.4 System turn off

Shut down MMI CELLTOOLS® by using the Exit command, in the File drop down menu.

File->Exit


Or press the cross on the right upper corner of the program window. Shut down the laser controller by pressing the button to deactivate the laser and then rotating the key to the off (vertical) position to deenergize the electronics controller.

Shutdown the PC as needed by using the keyboard command.

6.3 MMI CELLTOOLS® – Main Application

6.3.1 Main Window and plugins

Extended tool area:



Extended tool area

The MMI CELLTOOLS® provides a live video screen of the actual field of view of the microscope. To get a fast overview of the sample, scan functions are available in the overview panel on the right upper corner of the program window.

Live video and overview scan

On the upper function bar the program gives access to the user specific database. This database handles and stores all the parameters defined by user.

User specific database

Each plug-in(or in other words module) appears as a separate sheet on the right side of the program window. To switch from one plug-in to the other user only has to click on the appropriate sheet.


Switching between plugins

6.3.2 User specific database


All settings saved in the MMI CELLTOOLS® are unique for the Microsoft Windows® user actually logged in. MMI CELLTOOLS® fully supports the Microsoft Windows® user management. During program start the last settings saved by the active user are loaded.

Full support of Microsoft Windows® user management

In the slide selection box you find all samples you defined in the past. All documentation is saved under this name. To



change the database you press the **E** button and get the slide Editor:




+ Define new slide, all parameters from the active slide will be copied. You rename the slide by clicking on the name.

- Remove the active slide

Navigation buttons Navigate through the slides

In the SetUp selection box you find all setups you defined for the active slide. To change the database you press the “E” button and get the SetUp Editor. To use this Editor see Slide Editor described above. If you create a new slide all setups will automatically be copied from the active slide.



If you run different experimental settings, e.g., microdissection with bright field, microdissection with fluorescence of FITC or tweezers it is recommended to define one setup for each of these situations.

In the Objective selection box you find all objectives you defined for the active setup. To change the database you press the “E” button and get the Objective Editor. To use this editor follow the method of Slide Editor described above. If you create a new setup all objectives will automatically be copied from the active slide.



The distance between two points in the visible sample area depends of the magnification of the active objective. Therefore the program needs information on the chosen magnification for the table control. From the menu, choose the appropriate objective for this magnification.

Choose and confirm magnification

If the xy-stage is still not following the mouse movement the objective has to be recalibrated (See Objective calibration).

6.3.3 Video settings


In the standard setup the MMI CELLTOOLS® are using a specific frame grabber board. Only with this frame grabber board is the fluorescence mode fully supported. In this case the user has no video setting available.

In special cases the MMI CELLTOOLS® also support Microsoft Windows® vfw-video drivers. If the standard board is not installed, the user can select between different video drivers supported by Windows®. Depending on the video driver selected specific video setting can be available.

6.3.4 Camera configuration

For different illumination setups, like bright field, phase contrast or fluorescence, the ccd camera needs different parameters for integration time and gain to produce good images. The camera settings can be changed in the menu

Settings->Camera RC:



For the gain the automatic adjustment by the camera can be selected. In this case the camera measures the light intensity and attempts to choose appropriate settings. With the **gain limit** you can define the maximum gain, which is used by the camera in automatic mode. Note that a high gain induces a lot

Auto illumination

of more noise and the image quality may become poorer.

For standard (i.e., white) light choose an exposure time from 1/125 to 1/100000 according to light conditions, sample and chosen objective. In the FL mode the camera corrects for the fluctuation in the illumination (50 Hz Europe, 60 Hz North America). The video screen shows you immediately the effects of the camera setting. For fine tuning use the scrollbar **Gain**. These settings are used until you change it again or switch to another objective or setup.

Standard light

For fluorescence mode choose an exposure from 0.1 up to 8 according to light conditions, sample and chosen objective. Typically, fluorochromes will differ in emitted intensity. For fine tuning use the scrollbar **Gain**. This setting is used until you change it again.

Fluorescence

As standard the camera control is set on Remote. In this mode MMI CELLTOOLS can switch the exposure time between standard and fluorescence. Other camera setups are part of the installation of the system. If you want to change to camera settings, which are not accessible through the camera menu of mmi CellTools, choose Manual.

Camera control

The recommended (and installed) settings for the ccd camera are listed in the annex of this manual.


If you change the illumination devices of your microscope you should run the automatic White balance. If the button is pressed the white balance is reset. If the button is released color settings will be used.

White balance

To run the white balance you have to observe a white object (white sheet of paper) with the microscope. If the observed object is not white, the white balance will fail and the standard settings will continue to be used.

6.3.5 Fluorescence

To define settings for a fluorescence experiment, it is recommended to start with defining a new setup (see chapter 6.3.2.) For a set-up name use "Fluorescence dye1", for example.

Now you can set the ccd camera on fluorescence exposure time by clicking the *Fluorescence* button . The fluorescence camera settings and video settings will automatically be started.

Long term integration CCD camera


To avoid photobleaching of the sample by the lamp intensity, MMI CELLTOOLS® offers you a frozen image to draw the cutting path. To freeze the image, click the right mouse button with the cursor over the video screen. Choose *Freeze* from the context menu. The video will be frozen, so that you can close the fluorescence illumination shutter. All procedures still work in the

Frozen image

frozen image, but you will see no results. To check the results you have to switch back to LTE (long term integration) mode by using the context related menu (right mouse button). After going back one step you have to open the fluorescence illumination shutter and check the results.

6.3.6 XY stage motion

The movement of the stage is controlled by the MMI CELLTOOLS®. Three different modes allow fast scanning for overview or slow scanning for details.

Choose  or open the context menu by pressing the right mouse button with the cursor in the video window and choose **Move** with a left mouse click. By pressing the left mouse button the table follows directly the mouse movement.

With the keyboard cursor keys the table can be moved slowly and continuous in the four directions.

The arrow keys in the number block on the PC keyboard move the sample step by step with a 10% overlap from one viewing area to the next. This allows searching for a particular location in the sample without missing any part of the sample.





Three table movement modes

Direct mouse control



Keyboard control

*Keyboard control:
step by step
movement*

6.3.7 Slide overview for navigation

In the right lower corner of the MMI CELLTOOLS® you will find a panel representing the area of the sample defined by  limit1 and  limit2 (See "Handling a new slide"). After selection of the area of interest with the  button and scanning the slide by pressing the start scan button , this window shows you an overview over your sample (the "roadmap image"). The position of the detail of the large video image on the slide is indicated with a red blinking point.

The maximum area of the scan is limited by the computer memory and the magnification of the microscope. As standard, a scan is executed with the 4x objective. In this case you can define the left upper corner of the inner window of the membrane slide as limit1 and the right lower corner as limit2. To do so, follow these steps:

- move the xy-stage to the left upper corner of the membrane
- press the limit1 button  and observe the small, blinking, red rectangle in the overview.
- move the xy-stage to the right lower corner of the membrane
- press the limit2 button  and observe the small, blink-




Limit definition



ing, red rectangle in the overview.



For higher magnifying objectives the limits must be selected closer together to get a good visible overview and not overload the PC memory. If the scan area is too big the white/red-blinking rectangle will be change to big red dot.

Define the actual area of the scan quickly with the  area tool. After pressing the button you can select the area of interest in the overview window using the PC mouse. Only the area of interest will be scanned.



Area of interest

With a right click on the area of interest you can define the image quality of your scan. The highest possible quality is depending depends on the PC resources.


Begin the scan by pressing  the Start Scan button. With the  stop scan button you can always interrupt the scan. Also pressing the **ESC** key will interrupt the scan.



start scan




Stop scan

With  you can open the overview image in the standard Windows® image viewer. The time that Windows® needs to open the application is dependent on the image size and PC performance.




WINDOWS®
image viewer

With the  copy button you can copy the overview image into the Windows® clipboard. Following this step the image is available for all other standard Windows® applications.



copy to clipboard

With the  save button you can save the overview image on your hard disk or other drive.



save image
on hard disk

A white/red-blinking rectangle represents the position of the detail visible in the large video window. You can move this red frame with the left mouse button. The motorized xy-table moves automatically to the chosen detail. With this navigation method you always see the position on the slide.

You can also move to the position of interest by double clicking into the overview area with the left mouse button.

Navigation


6.3.8 Additional Settings for the xy-stage

*Deceleration of
the xy-stage*

If you press in the menu bar item

Settings->XYtable >XYmotionparam

it takes following window.



In this window you can limit the speed and acceleration of your motorized xy-stage. This could be necessary if you mount a heavy load on your xy-stage or if you're working with liquid cells otherwise no changes are recommended from the factory settings..

Additionally you can define the time for a single laser shot (only for MMI UVCUT®), by default is 50 msec. A range of 50 – 100 msec. is typical.

*Single Laser pulse
duration*

6.3.9 Calibration procedures

The MMI CELLTOOLS® are calibrated by MMI during the instrument installation. When you remove or add objectives or if the original camera position has been changed you have to repeat the calibration as described below.

6.3.9.1 Camera orientation

To check the camera orientation:

- Using the 40x objective move one small and remarkable point of the sample to the left upper corner of the video panel.
- Now press the right arrow key in the number pad.

If you move the object horizontally to the right upper corner and the object is still on the same level, the camera orientation is useable. If the result shows a shift between the two corners, loosen the hex screws on the camera port and rotate the camera carefully. Repeat this "move and observe" operation again until

the object moves on the same level across the monitor view. When the object moves in a straight line the camera orientation is correct. Retighten screws.

6.3.9.2 Objective calibration

When you add a new objective to your microscope you have to calibrate the xy-stage control. Over the time of operation it may be necessary to recalibrate the control, when contours and cutting lines are no longer exactly fitting.

Take a slide with sample. Choose an objective in the objective revolver and select the corresponding objective in the program. If the objective does not exist create it as described in Section 6.3.2.

Open in the menu bar item


UVCut->Calibration->Objective calibration

Calibration for objectives



Confirm objective

The objective calibration window appears.



To start the calibration of the active objective you can press the button **Calibrate**.

The Calibration window will be minimized and you will find two marked areas in the top left and bottom right of the screen (calibration areas):

Space key

- Move one small and remarkable point of the sample in the upper left calibration area
- Mark it with a left mouse button click
- Move this point to the bottom right calibration area

- Mark this point again with a left button click
- Confirm calibration with a right mouse click.

Now the calibration window comes back. Accept the calibration by pressing **OK** or discard your changes by pressing **Cancel**.

Remark: You always can switch between the move mode and the mark mode by pressing the **Space** key on your keyboard. In the move mode you see the **hand courser** in the mark mode the **marker cursor**.

*Lens OffSet
calibration***6.3.9.3 Lens OffSet calibration**

When you install a new objective into your microscope you have to calibrate Lens OffSet. Over the time of operation it may be necessary to recalibrate the control, when contours and cutting lines are no longer exactly fitting.

Calibration:

1. Take a slide with sample.
2. Make an objective calibration if necessary (See Objective calibration).
3. Choose an objective in the objective revolver and select the corresponding objective in the program. Make a shape around some noticeable object. Move a noticeable object into center of field of view.
4. Change an objective (physically and programmatically).

Please note, the Lens offset calibration procedure should go from objective with smaller magnification to objective with higher one for example from 4x to 10x, and should have no gap like 4x to 20x instead of 4x to 10x or 10x to 20x etc.

Press the menu item

Settings->Lens OffSet calibration

Right mouse click to set up a difference between a shape and corresponding object on the screen.

Repeat steps 1-4 for each objective.

6.3.10 Multi-User Report

Multi-User Report


User can get information about time spend for mmi Cell Tools for each Windows user account.

Press the menu item

Settings->MultiUser Report

Or start it from

Start->Progreಾಮms->mmiCellTools-> MultiUser Report




N	User Name	Spend time
1	LB8	0000:07:28
2	Stefan	0001:06:07
3	hardware	0002:22:11

a time format is hhhh:mm:ss

6.3.11 Hardware self monitoring

*Hardware
self monitoring*

When system has a problem with some hardware, there is a red blinker. The detailing of problem (or problems) see in that combo box:




6.3.12 Overview Image Quality

*Overview Image
Quality*


For some special cases, for example, when you scanning a huge image in overview window under 20x or more for full slide, you can get a follow message:



for overcome these difficulties right mouse click under overview window



and choose *Overview Image Quality item*. Because of smart and fast algorithm, you will see nearly no visibly differences for 10 or yet more percents, but size will be match less.




by default *Image Quality* is 50%.

6.4 MMI CELLTOOLS® – UV CUT® Plug-in


The MMI UV CUT® plug-in appears as a separate sheet on the right side of the program window.

The MMI UVCut® plug-in



To switch from one plug-in to the other you only have to click on the appropriate sheet.


6.4.1 Drawing the cutting contour





The laser cuts along the contour you have drawn around the interesting area you want to microdissect.

You can choose between freehand drawing and preset forms such as circles, ellipses, lines and rectangles respectively.

Drawing modes

To activate the freehand drawing mode press . The cursor appears as a pencil. Press the left mouse button to draw the contour around the area you are interested in. The program closes the contour automatically if you have selected the option "AutoEnclose Shape" in the menu **UVCut->Draw and Cut**.

Freehand drawing

Choose  or . The cursor appears as a symbol of the chosen form. Move the cursor (pressed left button) over the length of the area you want to mark. If you right click on these buttons you can choose to draw a circle or ellipse or respectively a rectangle or line. You can also give a defined dimension to your mark by typing the value in the corresponding input box.

Preset form contours


The drawing mode can be activated through the buttons in the MMI UV CUT® plug-in or in the context menu. The context menu can be activating by pressing the right mouse button (with the cursor in the video window) and choose **Paint**. The last used drawing mode will be activated.

Switching to drawing mode

A quick switching between moving and drawing mode is possible by pressing the **Space** key on the keyboard.

You can delete the current contour by pressing the key. If you want to delete any object choose the delete mode with

Delete contours

 button and click on the contour you want to delete.

With the context menu you also can delete one or all contours of a group or all contours on your slide.

6.4.2 Using the MMI IsolationCap®

The MMI IsolationCap® for a Single Step Collection contains a diffuser insert in the lid. Position the lid over the area of tissue you want to review for maximum image quality during microdissection operations.

To mount the tube simply slide the lid into the cap holder and open the tube. The tube holder can easily be fixed to the automatic cap lift. The lift has a small magnet and two positioning pins to hold the Inverted tube.


You can cut an object with the Caplift up or down. But the MMI IsolationCap® works most reliably if microdissection is done with the cap lowered onto the membrane surface.

The xy-stage can also be moved when the cap is down. The maximum moveable distance, without lifting the cap, should be limited to avoid tearing or other damage to the sample membrane. The Caplift free distance is 1000 microns (1mm) as the default. Under specific conditions it could be helpful to increase the free distance to be able to cut bigger objects with the cap down.

The Caplift free distance can be changed by using the menu item

UVCut->Caplift free distance

The Caplift free distance window appears:



If the cap is down on the tissue the cap will automatically be lifted whenever you leave the free distance area by moving the xy-stage. Then it will be repositioned automatically.

Image improvement using the mmi Isolation-Cap®



Caplift Up/Down button

Automatic cap up/down function

cally.

If the Caplift is down you can cut only objects fitting into the free distance area. If you work with a 0.5ml cap it is recommended to set a distance of about 1-2mm. The free distance can be adjusted as needed when working with the 0.2ml or 1.5ml tubes.

6.4.3 Cutting of objects

The automated cutting along the contours makes your work easy and fast. You can decide if you want to cut one area or one group of areas.

For most applications the best cutting results are observed with the 20x or the 40x objective. The user must select the appropriated objective in the objective selection box.

The laser set up includes the stage's parameters: **speed**, **focus** and **power**. It is essential for the cutting performance to properly set these settings for each objective and tissue type.

The **speed** slider defines how fast the xy-stage moves during the cutting procedure itself. It does not affect any other movement parameters of the xy-stage.

The **focus** slider adjusts the focus position of the laser beam in z-direction within the tissue sample. The cutting performance is very sensitive to this parameter. For example, the user would expect to use different focus settings for soft tissue vs. bone.

The **power** slider defines the laser power. The power needed at the sample normally will be proportional to sample thickness. Choose the power setting, which enables clean laser cutting. In laser microdissection, the laser is used as a "light scalpel". The correct combination of focus and power define the sharpness of the scalpel tip, i.e., the overall cutting line width for a given speed. These are the most important parameters for the user to master.

Before dissecting the tissue of interest, test the laser performance in another area of the sample or on a sacrificial slide with the same tissue. Draw a line or a circle and start cutting this figure with a low speed. Change the focus and power parameters during cutting to observe the effects by the laser. The sharpest line widths can be reached with power as low as possible and exact focus.


To adjust the parameters start with high power value. The focus can be adjusted easier with high power. Make several cuts with different focus values and compare the results. With a reduced power the focus will be in a more defined area. Decrease the power step by step and repeat focusing until the cutting line will be continuous, fine and clear.

Restriction for object size if cap lift is down

Automated cutting



Magnification



Laser set up

Test of cutting

Adjustment procedure


A cutting timer is selectable in the menu bar under

UVCut->Draw and Cut->Time left window

This timer gives you an estimate of the cutting time. To get access to the sliders click on the timer with the left mouse button. The timer will disappear and comes back if you do not use the sliders for more than three seconds.

Timer and access to cutting parameters during cutting

You can store as many of the laser parameter settings that you want for each objective. To define new setups see Chapter 6.3.2.



Laser parameter settings

Laser microdissection can only be started with .



Start cutting

„Repeat“ field does number of times of cutting.




Number of times


6.4.4 Group handling

If you have several types of cells to cut, you can define an unlimited number of groups. If you draw an object it will always be assigned to the actual group shown in the group selection box.

The most common reason for using more than one group is to permit the collection of different types of cells in different MMI IsolationCaps® (i.e., Group A with Cap A, Group B with Cap B, etc.).


All drawn areas belonging to one group are marked with the same contour color.

You can adjust the drawing attributes, such as group name, color and line thickness in the editor by pressing the  button in the group name box. The group editor window appears:




Several dissections in one operation

Contour color and thickness





6.4.5 Distance measurement

Click  to measure different distances on the screen. The cursor appears as a ruler. Press the left mouse button and drag the mouse. After releasing the mouse button you precisely see the measured distance.



Distance measurement tool

6.4.6 Pin positions

With the  button you can save the current xy-stage position. With the  button you can delete it again. You can save as many positions as you want.



create and delete pin positions


You can move the xy-stage from one pin position to the next by moving the pin position slider.

6.4.7 Autodocumentation

To set the auto documentation features, open the auto documentation window from the menu bar:

UVCut->Autodocumentation->Set Autodocumentation

Auto documentation for images before and after each cut



In the window you can activate the auto documentation and define the parameters of the function. The program creates automatically:

- Text Report (.htm file)
- First Image
- Image for each shape of a group (if not marked only an image before and after cutting a whole group will be saved)
- Last image

Note that you should work with a jpeg image in order to avoid accumulation of many, large TIFF or BMP files on the PC hard drive. The images in tiff format don't appear in the report.

6.5 Calibration

6.5.1 Laser position

Precise cutting requires that the laser beam position be exactly marked on the viewing screen. Over time and due largely to mild vibration, to temperature extremes in the laboratory environment, or to accidental moving or jarring of the instrument, recalibration of the laser position may be required.

Laser position

If the cutting action no longer follows the marked line you should adjust the laser position. It should be done for each magnification separately.



Cut a hole in the sample with . A red arrow will illuminate on the button when the laser fires or is on constantly.

Open **UVCut** from the menu bar and choose “**Set Laser Position**” in the menu “**Calibration**”. Fit the cursor cross with the hole, which had been cut and confirm this position with a right mouse button click. Continue with the next objective until the beam position has been defined for each one.

6.5.2 Correction of slide offset

If you have removed the slide from the xy-stage that you have drawn may not fit exactly when viewed on the monitor, if the slide is returned to the xy-stage for further work. [0]For correction, you may reposition the contours by pressing and holding the **Alt** key as you move the stage with the mouse. When the contours once again fit properly to the image, release the **Alt** key.

7 Disposables

7.1 Slides

Part No	Description
50101	RNase free, individually wrapped
50102	RNase free, individually wrapped, Box of 50 pieces
50103	Boxed, not individually wrapped, Box of 50 pieces
50104	Membrane, A4 sheet
50105	Glue Fixogum, 125 g

7.2 Caps

Part No	Description
50201	tube with adhesive lid, with diffusor, 500 µl (standard size)
50202	tube with adhesive lid, with diffusor, 500 µl (standard size) pack of 50 pieces
50203	tube with adhesive lid, without diffusor, 500 µl (standard size)
50204	tube with adhesive lid, without diffusor, 500 µl (standard size) pack of 50 pieces
50205	tube with adhesive lid, with diffusor 200 µl
50206	tube with adhesive lid, with diffusor 200 µl, pack of 50 pieces
50207	tube with adhesive lid, without diffusor 200 µl
50208	tube with adhesive lid, without diffusor 200 µl, pack of 50 pieces
50209	tube with adhesive lid, with diffusor 1.5 ml
50210	tube with adhesive lid, with diffusor 1.5 ml, pack of 50 pieces
50211	tube with adhesive lid, without diffusor 1.5 ml
50212	tube with adhesive lid, without diffusor 1.5 ml, pack of 50 pieces

7.3 Cell Chamber

Part No	Description
50301	cell chamber with membrane complete with Petri dish coated with silicone pack of 10 pieces

8 Maintenance

Maintenance and service should be only performed by qualified MMI personnel or our designated representative(s). The MMI CELLCUT® contains no user-serviceable parts.

*Maintenance and
service*

9 Trouble shooting

This Section provides help for problems that can occur when working with UV CUT. Most of the problems can be easily solved. When needed, please refer to the problem description and the trouble shooting tips.

9.1 View

Image not clear, too dark, too bright or no image on monitor

View

- Check the camera set up and compare it with the camera manual.
- Check that there is sufficient illumination from the white light, and that the light is not obstructed.
- For improved image quality, use the diffuser.
- The microscope port button should be set on “side”.

9.2 Movement

The movement is not accurate and repeatable; the stored positions cannot be found again.

Movement

- Check if the sample is fixed correctly on the stage and cannot move.

9.3 Drawing

The drawing line is difficult to see.

Drawing

- Choose dark line colors for bright samples and bright line colors for dark samples. Increase the thickness of the line if needed.

9.4 Cutting

No cutting can be observed.

Cutting

- The laser has not been switched on (key and button).
- The magnification chosen not correspond to the objective in use.
- The laser power is too low.
- The laser focus is not well adjusted.
- The cutting is not perfect; bridges remain in the cutting path. Repeat the cut if necessary.
- The cutting speed is too high.
- The cutting path does not correspond to the drawn contour. Check the camera orientation (6.3.9.1). The laser position

may need to be adjusted. Incorrect magnification may be selected.

- The calibration may have changed. Recalibrate when necessary.

9.5 Collection

The dissected objects are not collected.

Collection

- Make sure that the lid is placed directly over the treated area and lowered down onto the membrane. Ensure that it is not on the metal frame.
- Assure that the cutting action has been complete and successful. Adjust laser settings or repetitions to achieve full penetration of the beam through the membrane and tissue.

10 Warranty

In Europe and Asia:

The optical system and the control system are protected against defects in materials and workmanship for a period of one year from the date of installation. This warranty is limited to repair or replacement of system components which returned to MMI within the given one (1) year period. Unauthorized repairs, replacements and modification will cancel all warranties.

software
optical system

The microscope is not included in the warranty given by MMI **AG**. The microscope warranty depends on the manufacturer of the microscope, the country of location, and the responsible distributor. Consult Service at MMI AG for details.

microscope

In North America:

For the standard CellCut instrument configuration, the optical system, including the microscope, and the control system are protected against defects in materials and workmanship for a period of one year from the date of installation. Any failure resulting from unauthorized system modification, abuse, neglect, misuse, or acts of nature are not covered under this or any warranty.

Custom instrument configurations may be covered under the terms of the standard warranty via written acknowledgement by MMI Inc.

Unauthorized instrument repairs, replacements, or modifications cancel all warranties expressed or implied.

Consult Service at MMI Inc. for details.

11 Customer support

For questions about the MMI CELLCUT®

contact

(technical, consumables, warranties) please contact:

Europe and Asia:



service@molecular-machines.com



Molecular Machines & Industries AG (MMI),
Flughofstrasse 37, CH - 8152 Glattbrugg,
Switzerland



+41 - (0)1 - 8 09 10 10



+41 - (0)1 - 8 09 10 11

North America:



MMI Molecular Machines & Industries Inc.
North America
Don Armstrong, Ph.D.
P.O. Box 23991
Knoxville, TN 37933-1991

USA



+1 - (865) 988 7500



+1 - (865) 988-6666

MMI via the Internet:

www.molecular-machines.com

12 System data overview

Required minimum workspace:

The table for the microscope, laser, optical equipment, computer monitor and keyboard requires a minimum dimension of 1,20m x 0,90m

The computer should be positioned under or near the table.
The camera/computer connection cable is 2m long to ensure reliable data transfer.

Components:

Name	Description	Specifications
microscope	research microscope, Zeiss; Nikon (standard configuration); Olympus; Leica	(inverted) with fluorescence port
xy scanning table	scanning table with stepper motors	scanning area : (120x100)mm ² accuracy : < 1µm
laser system		wavelength : 355nm; line voltage : 200-240 VAC or 100-110 VAC, 50-60 Hz, 1,0A
Optical and electronic coupling unit	Laser control unit MMI	
ccd camera with power supply	DXC-390P	3CCD Color Video Camera; CCD-Chip : 1/3-inch 3-Chip with Power HAD technology horizontal resolution : 800 TV lines; Sensitivity : 2000 lux at F8; single to noise ratio : 61dB
computer	Pentium PC with Win XP Pro	With RGB frame grabber
control software	MMICELLTOOLS with MMIUVCUT® plug-in	32 bit Software

A Software and hot key List

A.1 Shortcuts:

MMI CELLTOOLS® application run **Ctrl+Alt+F12**

MMI CELLTOOLS® Main Window

Main Menu

File

	Alt+F
New Slide	Ctrl+N
Save Image	Ctrl+S
Start Record Movie	Ctrl+Alt+S
Stop Record Movie	Ctrl+Alt+D
Exit	Alt+F4

Edit

	Alt+E
Copy Image (for Move mode)	Ctrl+C
Copy Shape (for Paint mode)	Ctrl+C
Insert Shape (for Paint mode)	Ctrl+V

Video

Driver (non Flashpoint card)	Alt+V
	Ctrl+D
Video Type (non Flashpoint card)	Ctrl+T

Help

Help Topics	Alt+H
	F1
Info	Alt+I

MMI UVCUT® only

Cap lift up/down	F2
Select next marker	Ctrl +Tab
Select prev marker	Shift+Tab

A.2 Hotkeys:***For any window***

Accept changes & close window	ENTER key
Close current opened window	ESC key
Stop current procedure or cancel last command	ESC key

Main Window

Changing type of cursor	SPACEBAR key
Screen to first marker go	HOME key
Screen to last marker go	END key
Screen to next marker go	+ key
Screen to prev marker go	- key
Screen to next marker go	PAGE UP key
Screen to prev marker go	PAGE DOWN key
Screen 20% go up	UP ARROW key
Screen 20% go down	DOWN ARROW key
Screen 20% go left	LEFT ARROW key
Screen 20% go right	RIGHT ARROW key
Screen 20% go up	MOUSE WHEEL UP
Screen 20% go down	MOUSE WHEEL DOWN
Screen 90% go left	NUMPAD 6 key
Screen 90% go right	NUMPAD 4 key
Screen 90% go up	NUMPAD 8 key
Screen 90% go down	NUMPAD 2 key
Emergency stop motion	NUMPAD 5 key

Shift the contours relative to the sample
(Press ALT while moving the sample)

ALT+ Move mode